

Fortify or Calibrate for 203 Pesticides by GC-MS/MS With Our New Certified Reference Materials Kit



- Accurately identify and quantify pesticide residues by GC-MS/MS in fruits, vegetables, botanicals, and herbals like tea, ginseng, ginger, Echinacea, and dietary supplements.
- Comprehensive 203-compound kit covers food safety lists by the FDA, USDA, and other global governmental agencies; individual ampuls also sold separately.
- Formulated and grouped for maximum long-term stability.
- Quantitatively tested to confirm composition; detailed support documentation provided.
- Optimized multiresidue pesticide method is offered free of charge; downloadable XLS file includes conditions and transition tables.
- Certified reference material (CRM) manufactured and QC-tested in Restek's ISO-accredited labs satisfies your ISO requirements.
- Restek is your complete supplier for world-class GC-MS/MS multiresidue pesticide analysis.

No more long nights or weekends in the lab. No more custom standards. Restek's food safety experts can help you make quick work of getting the accurate results you need. Combine this ready-made multiresidue pesticide standards kit with Restek's internal standards, Rxi®-5ms GC columns, Q-sep® QuEChERS sample preparation, Sky® inlet liners, accessories, and more.

Learn more at www.restek.com/GC-multiresidue

Perform Quantitative GC-MS/MS Multiresidue Pesticide Analysis of 203 Compounds With 1 Reference Standards Kit (cat.# 32562)

Std #**	Compound (100 µg/mL each)	CAS #	Std #**	Compound (100 µg/mL each)	CAS #	Std #**	Compound (100 µg/mL each)	CAS #
3	2,3,5,6-Tetrachloroaniline	3481-20-7	2	Dieldrin	60-57-1	4	Oxadiazon	19666-30-9
2	2,4'-DDD	53-19-0	4	Dimethachlor	50563-36-5	3	Oxyfluorfen	42874-03-3
2	2,4'-DDE	3424-82-6	4	Diphenamid	957-51-7	5	Paclobotrazol	76738-62-0
2	2,4'-DDT	789-02-6	3	Diphenylamine	122-39-4	9	Parathion (Ethyl parathion)	56-38-2
2	2,4'-Methoxychlor	30667-99-3	9	Disulfoton	298-04-4	4	Pebulate	1114-71-2
3	2,6-Dichlorobenzonitrile (Dichlobenil)	1194-65-6	8	Edifenphos	17109-49-8	5	Penconazole	66246-88-6
7	2-Phenylphenol	90-43-7	2	Endosulfan ether	3369-52-6	3	Pendimethalin	40487-42-1
3	3,4-Dichloroaniline	95-76-1	2	Endosulfan I	959-98-8	3	Pentachloroaniline	527-20-8
2	4,4'-DDD	72-54-8	2	Endosulfan II	33213-65-9	2	Pentachloroanisole	1825-21-4
2	4,4'-DDE	72-55-9	2	Endosulfan sulfate	1031-07-8	2	Pentachlorobenzene	608-93-5
2	4,4'-DDT	50-29-3	2	Endrin	72-20-8	3	Pentachlorobenzonitrile	20925-85-3
2	4,4'-Dichlorobenzophenone	90-98-2	2	Endrin aldehyde	7421-93-4	3	Pentachloronitrobenzene (Quintozene)	82-68-8
2	4,4'-Methoxychlor olefin	2132-70-9	2	Endrin ketone	53494-70-5	2	Pentachlorothioanisole	1825-19-0
7	Acequinocyl	57960-19-7	1	EPN	2104-64-5	6	Phenothrin (<i>cis</i> & <i>trans</i>)	26002-80-2
4	Acetochlor	34256-82-1	3	Ethalfuralin	55283-68-6	9	Phorate	298-02-2
6	Acrinathrin	101007-06-1	8	Ethion	563-12-2	1	Phosalone	2310-17-0
4	Alachlor	15972-60-8	2	Ethylan (Perthane)	72-56-0	1	Phosmet	732-11-6
2	Aldrin	309-00-2	5	Etofenprox	80844-07-1	9	Piperonyl butoxide	51-03-6
4	Allidochlor	93-71-0	5	Etridiazole	2593-15-9	1	Pirimiphos ethyl	23505-41-1
2	alpha-BHC	319-84-6	8	Fenamiphos	22224-92-6	1	Pirimiphos methyl	29232-93-7
6	Anthraquinone	84-65-1	5	Fenarimol	60168-88-9	4	Pretilachlor	51218-49-6
5	Atrazine	1912-24-9	8	Fenchlorphos (Ronnel)	299-84-3	4	Prochloraz	67747-09-5
1	Azinphos ethyl	2642-71-9	1	Fenitrothion	122-14-5	5	Procymidone	32809-16-8
1	Azinphos methyl	86-50-0	4	Fenpropathrin	39515-41-8	3	Prodiamine	29091-21-2
3	Benfluralin	1861-40-1	2	Fenson	80-38-6	8	Profenofos	41198-08-7
2	beta-BHC	319-85-7	8	Fenthion	55-38-9	3	Profluralin	26399-36-0
6	Bifenthrin	82657-04-3	6	Fenvalerate	51630-58-1	4	Propachlor	1918-16-7
6	Bioallethrin	584-79-2	5	Fipronil	120068-37-3	4	Propanil	709-98-8
3	Biphenyl	92-52-4	7	Fluzifop- <i>p</i> -butyl	79241-46-6	5	Propargite	2312-35-8
8	Bromfeninfos-methyl	13104-21-7	3	Fluchloralin	33245-39-5	4	Propisochlor	86763-47-5
8	Bromfenvinphos	33399-00-7	6	Flucythrinate	70124-77-5	4	Propyzamide	23950-58-5
8	Bromophos ethyl	4824-78-6	5	Fludioxonil	131341-86-1	8	Prothiofos	34643-46-4
8	Bromophos methyl	2104-96-3	4	Fluquinconazole	136426-54-5	1	Pyraclufos	77458-01-6
7	Bromopropylate	18181-80-1	5	Fluridone (Sonar)	59756-60-4	1	Pyrazophos	13457-18-6
5	Bupirimate	41483-43-6	5	Flusilazole	85509-19-9	4	Pyridaben	96489-71-3
5	Captafol	611/2425	4	Flutolanil	66332-96-5	1	Pyridaphenthion	119-12-0
5	Captan	133-06-2	5	Flutriafol	76674-21-0	5	Pyrimethanil	53112-28-0
8	Carbophenothion	786-19-6	5	Folpet	133-07-3	5	Pyriproxyfen	95737-68-1
7	Carfentrazone ethyl	128639-02-1	9	Fonofos	944-22-9	1	Quinalphos	13593-03-8
2	Chlorbenside	103-17-3	2	gamma-BHC (Lindane)	58-89-9	6	Resmethrin	10453-86-8
5	Chlorfenapyr	122453-73-0	2	Heptachlor	76-44-8	8	Sulfotepp	3689-24-5
2	Chlorfenson (Ovex)	80-33-1	2	Heptachlor epoxide (Isomer B)	1024-57-3	8	Sulprofos	35400-43-2
8	Chlorfenvinphos	470-90-6	2	Hexachlorobenzene	118-74-1	6	tau-Fluvalinate	102851-06-9
7	Chlorobenzilate	510-15-6	5	Hexazinone (Velpar)	51235-04-2	5	Tebuconazole	107534-96-3
2	Chloroneb	2675-77-6	8	Iodofenphos	18181-70-9	4	Tebufenpyrad	119168-77-3
3	Chlorothalonil	1897-45-6	5	Iprodione	36734-19-7	6	Tefluthrin	79538-32-2
7	Chlorpropham	101-21-3	1	Isazophos	42509-80-8	5	Terbacil	5902-51-2
1	Chlorpyrifos	2921-88-2	2	Isodrin	465-73-6	8	Terbufos	13071-79-9
1	Chlorpyrifos methyl	5598-13-0	3	Isopropalin	33820-53-0	5	Terbuthylazine	5915-41-3
8	Chlorthiophos	60238-56-4	6	lambda-Cyhalothrin	91465-08-6	3	Tetrachloronitrobenzene (Tecnazene)	117-18-0
7	Chlozolinate	84332-86-5	5	Lenacil	8/1/2164	8	Tetrachlorvinfos	22248-79-9
2	<i>cis</i> -Chlordane	5103-71-9	8	Leptophos	21609-90-5	2	Tetradifon	116-29-0
2	<i>cis</i> -Nonachlor	5103-73-1	4	Linuron	330-55-2	6	Tetramethrin	7696-12-0
6	<i>cis</i> -Permethrin	61949-76-6	8	Malathion	121-75-5	3	THPI (Tetrahydrophthalimide)	1469-48-3
4	Clomazone (Command)	81777-89-1	7	Metalaxyl	57837-19-1	8	Tolclofos-methyl	57018-04-9
8	Coumaphos	56-72-4	4	Metazachlor	67129-08-2	3	Tolyfluanid	731-27-1
4	Cycloate	1134-23-2	8	Methacrifos	62610-77-9	2	<i>trans</i> -Chlordane	5103-74-2
6	Cyfluthrin	68359-37-5	4	Methoxychlor	72-43-5	6	Transfluthrin	118712-89-3
6	Cypermethrin	52315-07-8	9	Methyl parathion	298-00-0	2	<i>trans</i> -Nonachlor	39765-80-5
5	Cyprodinil	121552-61-2	4	Metolachlor	51218-45-2	6	<i>trans</i> -Permethrin	61949-77-7
7	DCPA methyl ester (Chlorthal-dimethyl)	1861-32-1	9	Mevinphos	7786-34-7	5	Triadimefon	43121-43-3
2	delta-BHC	319-86-8	5	MGK-264	113-48-4	5	Triadimenol	55219-65-3
6	Deltamethrin	52918-63-5	2	Mirex	2385-85-5	4	Triallate	2303-17-5
4	Diallate (<i>cis</i> & <i>trans</i>)	2303-16-4	5	Myclobutanil	88671-89-0	9	Triazophos	24017-47-8
1	Diazinon	333-41-5	4	N-(2,4-Dimethylphenyl)formamide	60397-77-5	5	Tricyclazole (Beam)	41814-78-2
3	Dichlofluanid	1085-98-9	3	Nitralin	4726-14-1	5	Triflumizole	68694-11-1
3	Dichloran	99-30-9	3	Nitrofen	1836-75-5	3	Trifluralin	1582-09-8
			4	Norflurazon	27314-13-2	5	Vinclozolin	50471-44-8

Figure 1: GC Multiresidue Pesticide Standard #1-OPP on Rxi®-5ms by GC-MS




Figure 2: GC Multiresidue Pesticide Standard #2-OCP on Rxi®-5ms by GC-MS




Figure 3: GC Multiresidue Pesticide Standard #3-ONP on Rxi®-5ms by GC-MS





Figure 4: GC Multiresidue Pesticide Standard #4-ONP on Rxi®-5ms by GC-MS



For notes, see page 5.

Figure 5: GC Multiresidue Pesticide Standard #5-ONP on Rxi®-5ms by GC-MS




Figure 6: GC Multiresidue Pesticide Standard #6-SPP on Rxi®-5ms by GC-MS




Figure 7: GC Multiresidue Pesticide Standard #7-HME on Rxi®-5ms by GC-MS




Figure 8: GC Multiresidue Pesticide Standard #8-OPP on Rxi®-5ms by GC-MS





Figure 9: GC Multiresidue Pesticide Standard #9-OPP on Rxi®-5ms by GC-MS



NOTES:

* For compound list, see page 2. For conditions, retention times, and transitions, visit www.restek.com/GC-multiresidue

** When combining a large number of compounds with different chemical functionalities, mix stability can be an issue. In formulating these standards, we extensively studied the 203 compounds involved, then grouped them into as few mixes as possible while still ensuring maximum long-term stability and reliability. For quantitative analysis, we recommend analyzing each mix separately to ensure accurate results for every compound.

For Multiresidue Pesticide Analysis by GC-MS/MS, Restek Recommends...

Certified Reference Materials (CRMs)

GC Multiresidue Pesticide Kit

- For a complete compound list, see page 2.

Contains 1 mL each of these mixtures.
cat.# 32562 (kit)



Restek is your complete supplier for world-class GC-MS/MS multiresidue pesticide analysis: reference and internal standards, Rxi®-5ms GC columns, Q-sep® QuEChERS sample preparation, Sky® inlet liners, accessories, and more!

GC Multiresidue Pesticide Standard #1 (16 components)

Organophosphorus Compounds
100 µg/mL each in toluene, 1 mL/ampul
cat.# 32563 (ea.)

GC Multiresidue Pesticide Standard #2 (40 components)

Organochlorine Compounds
100 µg/mL each in toluene, 1 mL/ampul
cat.# 32564 (ea.)

GC Multiresidue Pesticide Standard #3 (25 components)

Organonitrogen Compounds
100 µg/mL each in toluene:acetonitrile (99:1), 1 mL/ampul
cat.# 32565 (ea.)

GC Multiresidue Pesticide Standard #4 (28 components)

Organonitrogen Compounds
100 µg/mL each in toluene, 1 mL/ampul
cat.# 32566 (ea.)

GC Multiresidue Pesticide Standard #5 (34 components)

Organonitrogen Compounds
100 µg/mL each in toluene, 1 mL/ampul
cat.# 32567 (ea.)

GC Multiresidue Pesticide Standard #6 (18 components)

Synthetic Pyrethroid Compounds
100 µg/mL each in toluene, 1 mL/ampul
cat.# 32568 (ea.)

GC Multiresidue Pesticide Standard #7 (10 components)

Herbicide Methyl Esters
100 µg/mL each in toluene, 1 mL/ampul
cat.# 32569 (ea.)

GC Multiresidue Pesticide Standard #8 (24 components)

Organophosphorus Compounds
100 µg/mL each in toluene, 1 mL/ampul
cat.# 32570 (ea.)

GC Multiresidue Pesticide Standard #9 (8 components)

Organophosphorus Compounds
100 µg/mL each in toluene, 1 mL/ampul
cat.# 32571 (ea.)

Labeled Pesticide Residue Internal Standards for Food Analysis

- Isotopically labeled to provide the best approach for pesticide residue quantification.
- Multiple options let you choose internal and surrogate standards that will mitigate matrix effects.
- Economically priced and compatible with both LC-MS and GC-MS applications; even helpful for optimizing LC-MS/MS system performance.
- Certified reference material (CRM) manufactured and QC-tested in Restek's ISO-accredited labs—satisfy your ISO requirements.

Compound	CAS #	Solvent	Conc.	cat.#
Atrazine-d5	163165-75-1	ACN	100	31984
Carbaryl-d7	362049-56-7	ACN	100	31985
Diazinon-d10 (diethyl-d10)	100155-47-3	ACN	100	31986
Dichlorvos-d6	203645-53-8	A	100	31987
Dimethoate-d6	1219794-81-6	ACN	100	31988
Diuron-d6	1007536-67-5	ACN	100	31989
Linuron-d6	330-55-2	ACN	100	31990

A = acetone; ACN = acetonitrile
Volume is 1 mL/ampul. Concentration is µg/mL.

QuEChERS Standards

- Ready to use for QuEChERS extractions—no dilutions necessary.
- Support for GC and LC with MS, MS/MS, and selective detectors.

Compatible with all major methods, including mini-multiresidue, AOAC, and European procedures. Save time with convenient mixes or make your own blend using our full line of single-component solutions.

QuEChERS Internal Standard Mix for GC-MS Analysis

(6 components)

PCB 18 (37680-65-2)	50 µg/mL
PCB 28 (7012-37-5)	50
PCB 52 (35693-99-3)	50
Triphenylmethane (519-73-3)	10
Triphenylphosphate (115-86-6)	20
Tris-(1,3-dichloroisopropyl)phosphate (13674-87-8)	50
In acetonitrile, 5 mL/ampul	
cat.# 33267 (ea.)	

Visit www.restek.com/quenchers for more QuEChERS standards.

QuEChERS Products

Restek Q-sep® QuEChERS Products

- Ready-to-use tubes, no glassware required.
- Preweighed, ultra-pure sorbents.
- Support original unbuffered, AOAC (2007.01), European (EN 15662), and mini-multiresidue QuEChERS methods.

QuEChERS methods are fast, easy, and cost-effective, and Restek Q-sep® products make QuEChERS procedures even simpler. All extraction salts, sorbents, and sample tubes are included—no specialized equipment or glassware is required. Prepare samples more efficiently with a complete line of QuEChERS supplies from Restek.



Q-sep® QuEChERS Extraction Salts

Fast, Simple Sample Prep for Multiresidue Pesticide Analysis

- Salt packets eliminate the need for a second empty tube to transfer salts.
- Go green by using packets with reusable tubes.
- Convenient and easy to use.

Description	Material	Methods	qty.	cat.#
Q-sep QuEChERS Extraction Kit (Original)	4 g MgSO ₄ , 1 g NaCl with 50 mL Centrifuge Tube	original unbuffered	50 packets & 50 tubes	23991
Q-sep QuEChERS Extraction Salt Packets Only (Original)	4 g MgSO ₄ , 1 g NaCl	original unbuffered	50 packets	23992
Q-sep QuEChERS Extraction Kit (EN)	4 g MgSO ₄ , 1 g NaCl, 1 g TSCD, 0.5 g DHS with 50 mL Centrifuge Tube	European EN 15662	50 packets & 50 tubes	26235
Q-sep QuEChERS Extraction Salt Packets Only (EN)	4 g MgSO ₄ , 1 g NaCl, 1 g TSCD, 0.5 g DHS	European EN 15662	50 packets	26236
Q-sep QuEChERS Extraction Kit (AOAC)	6 g MgSO ₄ , 1.5 g NaOAc with 50 mL Centrifuge Tube	AOAC 2007.01	50 packets & 50 tubes	26237
Q-sep QuEChERS Extraction Salt Packets Only (AOAC)	6 g MgSO ₄ , 1.5 g NaOAc	AOAC 2007.01	50 packets	26238
Empty 50 mL Centrifuge Tube, Polypropylene			50-pk.	26239
Empty 50 mL Centrifuge Tube, FEP			2-pk.	23997

TSCD—trisodium citrate dihydrate

DHS—disodium hydrogen citrate sesquihydrate

NaOAc—sodium acetate

Q-sep® QuEChERS dSPE Tubes for Extract Cleanup

Fast, Simple Sample Prep for Multiresidue Pesticide Analysis

- Packaged in foil subpacks of 10 for enhanced protection and storage stability.

Multiple sorbents are used to extract different types of interferences.

MgSO ₄	removes excess water
PSA	removes sugars, fatty acids, organic acids, and anthocyanine pigments
C18	removes nonpolar interferences
GCB	removes pigments, sterols, and nonpolar interferences

Description	Methods	qty.	cat.#
2 mL Micro-Centrifuge Tubes Prefilled with dSPE Materials for Cleanup (1 mL Extract)			
150 mg MgSO ₄ , 25 mg PSA	original unbuffered, mini-multiresidue, European EN 15662	100-pk.	26215
150 mg MgSO ₄ , 25 mg PSA, 25 mg C18	mini-multiresidue	100-pk.	26216
150 mg MgSO ₄ , 25 mg PSA, 2.5 mg GCB	mini-multiresidue, European EN 15662	100-pk.	26217
150 mg MgSO ₄ , 25 mg PSA, 7.5 mg GCB	mini-multiresidue, European EN 15662	100-pk.	26218
150 mg MgSO ₄ , 50 mg PSA	AOAC 2007.01	100-pk.	26124
150 mg MgSO ₄ , 50 mg PSA, 50 mg C18	AOAC 2007.01	100-pk.	26125
150 mg MgSO ₄ , 50 mg PSA, 50 mg GCB	AOAC 2007.01	100-pk.	26123
150 mg MgSO ₄ , 50 mg PSA, 50 mg C18, 50 mg GCB	AOAC 2007.01	100-pk.	26219
150 mg MgSO ₄ , 50 mg C18	NA	100-pk.	26242
150 mg MgSO ₄ , 50 mg PSA, 50 mg C18, 7.5 mg GCB	universal	100-pk.	26243
15 mL Centrifuge Tubes Prefilled with dSPE Materials for Cleanup (6 mL and 8 mL Extract)			
1,200 mg MgSO ₄ , 400 mg PSA	AOAC 2007.01	50-pk.	26220
1,200 mg MgSO ₄ , 400 mg PSA, 400 mg C18	AOAC 2007.01	50-pk.	26221
1,200 mg MgSO ₄ , 400 mg PSA, 400 mg C18, 400 mg GCB	AOAC 2007.01	50-pk.	26222
1,200 mg MgSO ₄ , 400 mg C18	similar to AOAC 2007.01	50-pk.	26244
900 mg MgSO ₄ , 150 mg PSA	original unbuffered, European EN 15662	50-pk.	26223
900 mg MgSO ₄ , 150 mg PSA, 15 mg GCB	European EN 15662	50-pk.	26224
900 mg MgSO ₄ , 150 mg PSA, 45 mg GCB	European EN 15662	50-pk.	26225
900 mg MgSO ₄ , 150 mg PSA, 150 mg C18	similar to European EN 15662	50-pk.	26226
900 mg MgSO ₄ , 300 mg PSA, 300 mg C18, 45 mg GCB	similar to European EN 15662	50-pk.	26245
900 mg MgSO ₄ , 300 mg PSA, 150 mg GCB	NA	50-pk.	26126

PSA—primary and secondary amine

GCB—graphitized carbon black

also available:



Resprep® Pesticide Residue Cleanup SPE Cartridges

www.restek.com/spe



Thomson SINGLE StEP® Filter Vials

www.restek.com/singlestep



Syringe Filters

www.restek.com/syringe-filters

Rxi®-5ms Columns (fused silica)


(low-polarity phase; Crossbond® diphenyl dimethyl polysiloxane)

- General-purpose columns for semivolatiles, phenols, amines, residual solvents, drugs of abuse, pesticides, PCB congeners (e.g., Aroclor mixes), solvent impurities.
- Most inert column on the market.
- Tested and guaranteed for ultra-low bleed; improved signal-to-noise ratio for better sensitivity and mass spectral integrity.
- Temperature range: -60 °C to 330/350 °C.
- Equivalent to USP G27 and G36 phases.

Description	temp. limits	qty.	cat.#
15 m, 0.25 mm ID, 0.25 µm	-60 to 330/350 °C	ea.	13420
30 m, 0.25 mm ID, 0.25 µm	-60 to 330/350 °C	ea.	13423



Patented



Dual Vespel® Ring Inlet Seals Washerless, Leak-Tight Seals

- Does not require a separate washer.
- Requires less torque to seal.
- Does not require retightening of reducing nut after several oven cycles.
- Extends column lifetime by preventing oxygen from reaching the column.
- Same price as the regular inlet seals with washers.
- Siltek® treatment provides enhanced inertness versus stainless steel.

Use a stainless steel seal for analyses of unreactive compounds. To reduce breakdown and adsorption of active compounds, use a gold-plated or Siltek®-treated seal.

for Agilent GCs

0.8 mm ID Dual Vespel Ring Inlet Seal	2-pk. cat.#	10-pk. cat.#	50-pk. cat.#
Gold-Plated	21240	21241	23418
Siltek-Treated	21242	21243	23419
Stainless Steel	21238	21239	23420

1.2 mm ID Dual Vespel Ring Inlet Seal	2-pk. cat.#	10-pk. cat.#
Gold-Plated	21246	21247
Siltek-Treated	21248	21249
Stainless Steel	21244	21245

for Thermo 1300/1310 GCs

0.8 mm ID Dual Vespel Ring Inlet Seal	2-pk. cat.#	10-pk. cat.#
Gold-Plated	22243	22244
Siltek-Treated	22247	22248

1.2 mm ID Dual Vespel Ring Inlet Seal	2-pk. cat.#	10-pk. cat.#
Gold-Plated	22245	22246
Siltek-Treated	22249	22250



also available:
Merlin Microseal Septa
www.restek.com/GC-Acc

RESTEK
Pure Chromatography

Questions about this or any other Restek® product?
Contact us or your local Restek® representative (www.restek.com/contact-us).

Restek® patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks in Restek® literature or on its website are the property of their respective owners. Restek® registered trademarks are registered in the U.S. and may also be registered in other countries.

© 2015 Restek Corporation. All rights reserved. Printed in the U.S.A.

www.restek.com



Lit. Cat.# FFTS2199-UNV

RESTEK

Sky® 4.0 mm ID Precision® Inlet Liner w/ Wool

For Agilent GCs equipped with split/splitless inlets

ID x OD x L	qty.	cat.#
Precision, Sky Technology, Borosilicate Glass with Quartz Wool		
4.0 mm x 6.3 mm x 78.5 mm	5-pk.	23305.5

Sky® 4.0 mm ID Single Taper Inlet Liner w/ Wool

For Agilent GCs equipped with split/splitless inlets

ID x OD x L	qty.	cat.#
Single Taper, Sky Technology, Borosilicate Glass with Quartz Wool		
4.0 mm x 6.5 mm x 78.5 mm	5-pk.	23303.5

Sky® 4.0 mm ID Precision® Inlet Liner w/ Wool

For Bruker/Varian GCs equipped with 1177 inlets

ID x OD x L	qty.	cat.#
Precision, Sky Technology, Borosilicate Glass with Quartz Wool		
4.0 mm x 6.3 mm x 78.5 mm	5-pk.	23328.5

Sky® 4.0 mm ID Single Taper Inlet Liner w/ Wool

For Bruker/Varian GCs equipped with 1177 inlets

ID x OD x L	qty.	cat.#
Single Taper, Sky Technology, Borosilicate Glass with Quartz Wool		
4.0 mm x 6.5 mm x 78.5 mm	5-pk.	23332.5

Sky® 4.0 mm ID Precision® Inlet Liner w/ Wool

For PerkinElmer Auto SYS™ and Clarus GCs equipped with split/splitless inlets

ID x OD x L	qty.	cat.#
Precision, Sky Technology, Borosilicate Glass with Quartz Wool		
4.0 mm x 6.2 mm x 92.1 mm	5-pk.	23450.5

Sky® 2.0 mm ID Inlet Liner w/ Wool

For PerkinElmer Auto SYS™ and Clarus GCs equipped with split/splitless inlets

ID x OD x L	qty.	cat.#
Straight, Sky Technology, Borosilicate Glass with Quartz Wool		
2.0 mm x 6.2 mm x 92.1 mm	5-pk.	23451.5

Sky® 4.0 mm ID Single Taper Inlet Liner w/ Wool

For Thermo TRACE 1300, 1310 GCs equipped with SSL inlets

ID x OD x L	qty.	cat.#
Single Taper, Sky Technology, Borosilicate Glass with Quartz Wool		
4.0 mm x 6.5 mm x 78.5 mm	5-pk.	23447.5

Sky® 3.5 mm ID Precision® Inlet Liner w/ Wool

For Shimadzu 17A, 2010, and 2014 GCs equipped with split/splitless inlets

ID x OD x L	qty.	cat.#
Precision, Sky Technology, Borosilicate Glass with Quartz Wool		
3.5 mm x 5.0 mm x 95 mm	5-pk.	23320.5

Sky® 3.5 mm ID Single Taper Inlet Liner w/ Wool

For Shimadzu 17A, 2010, and 2014 GCs equipped with split/splitless inlets

ID x OD x L	qty.	cat.#
Single Taper, Sky Technology, Borosilicate Glass with Quartz Wool		
3.5 mm x 5.0 mm x 95 mm	5-pk.	23336.5

* 100% SATISFACTION GUARANTEE: If your Sky® inlet liner does not perform to your expectations for any reason, simply contact Restek® Technical Service or your local Restek® representative and provide a sample chromatogram showing the problem. If our GC experts are not able to quickly and completely resolve the issue to your satisfaction, you will be given an account credit or replacement product (same cat.#) along with instructions for returning any unopened product. (Do not return product prior to receiving authorization.) For additional details about Restek's return policy, visit www.restek.com/warranty