

# Drying Methods and Processes

## Samples Types and Preparation

Size of the sample and its preparation shall provide the following: sample structure homogeneity, short time of drying, good repeatability of measurement within a measurement series and drying process result comparable to a reference result (standardised method).



**Solid bodies**  
- light sample,  
- requires grinding



**Samples of low humidity**  
- heavy sample  
(10 - 15 g)



**Liquids**  
- increasing active surface of evaporation is recommended

## Sample Structure Transformation

Unfavorable physical processes, affecting the sample structure, may occur during the drying process. These are caused by dispersion of indications and mistakes made when undertaking assessment of actual moisture content for a particular sample.



**Crust Formation**  
It is a process where an impermeable layer is formed on a sample surface. This makes removal of humidity from the sample impossible. As a result the indication being an outcome of an analysis is lower than the sample reference value.



**Sample Burning**  
Such a process is a consequence of too high drying temperature. It results in a change of sample colour. When sample burning occurs then the sample humidity value is greater than its reference value.



**Heat Absorption**  
Dark in colour samples absorb more heat than the light ones. This accounts for application of lower drying temperatures while drying light in colour samples. Tests need to be carried out in order to select the right temperature value.

## Automatic Control of Sample Weight

Obtaining optimal results for a drying process depends on samples quantity and weight. Too heavy sample lengthens duration of the drying process. Too light sample works against repeatability of the results. This proves that control of the sample weight is inevitable.



- A WEIGHT CONTROL** (checkweighing thresholds)
- B GRAPH** (drying process curve, registered for dynamic state)
- C DRYING PROGRAMS** (drying mode type, temperature, auto switch-off)
- D SAMPLES** (name, code, end mass, tolerance, etc.)
- E SAMPLE DESCRIPTION** (information on how to prepare a sample for a drying process)

Initial mass	2.7548 g
0:00:10	0.1503% M
0:00:20	0.6258% M
----- Auto switch off -----	
0:08:08	Auto 1
Result	11.4789% M
----- Auto switch off -----	
0:11:05	Auto 2
Result	11.9058% M
----- Auto switch off -----	
0:13:55	Auto 3
Result	12.0502% M
----- Auto switch off -----	
0:15:20	Auto 4
Result	12.0858% M

Demonstrative printout of TEST function.

## Methods of Analysis Completion

TEST function analyses weight variation for a particular sample, occurring during the drying process. There are 5 different options for automatic shutdown. The user shall select option allowing him to end the drying process in a way ensuring that the obtained humidity value is as close as possible to a reference value. Among other automatic shutdown possibilities there are time-defined, manual and user-defined options, the latest one being the best adjusted to a particular sample characteristics.

Obtained humidity value depends on start mass of a particular sample – an optimal mass value shall be selected prior running the tests.

Respective method for a particular sample drying, shall be selected based on tests optimizing the sample size, drying temperature and method of analysis end procedure.



# Reports and Statistics

## Report on Drying Process

RADWAG moisture analyzers allow the user to make self-configured reports. Analysis, summaries, etc. may be printed by means of any office printer (PCL).


The report comprises three sections: the header (A), the data area (B) and the footer (C). Each section can be freely configured by a user.

### Drying Process Statistics

Moisture content analyses performed for the same sample are used to determine the sample moisture content variation within a specified time interval (Trend).

Example of a simple drying report, generated by MA R moisture analyzer.

Drying	
Date	05.01.2021
Time	6:32:18
Operator	Admin
Product	Prod-01
Program	MAR-1
Drying profile	Standard
Drying profile parameters	90 °C
Finish mode	Manual
Start mass	0.674 g
-----	
0:00:30	
0:01:00	
0:01:30	
0:02:00	
0:02:30	
0:03:00	
-----	
Status	Completed
Drying time	0:03:00
End mass	0.499 g
Result	25.964% M



Drying	
Start date	2021.01.28
Start time	11:34:44
Operator	Admin
Product	Corn
Drying mode	Corn-PRG01
Drying mode	Standard
Drying mode parameters	100 °C
Auto switch-off mode	Auto
Finish mode parameters	1 mg/60 s
Printout interval	0:00:30
Start mass	0.590 g
-----	
Date and time	2021.01.28 11:35:14
Drying time	0:00:30
Product	Corn
Current result	5.085% M
0:00:30	5.085% M
Humidity content	5.085% M
Dry mass content	94.915% D
Humid / Dry	5.357% R
Tare	0.007 g
Gross	0.567 g
Set temperature	100 °C
Current temperature	99 °C
-----	
Date and time	2021.01.28 11:35:44
Drying time	0:01:00
Current result	7.795% M
-----	
Date and time	2021.01.28 11:39:14
Drying time	0:04:30
Current result	14.237% M
-----	
Status	Completed
End date	2021.01.28
End time	11:39:22
Drying time	0:04:38
Operator	Admin
Product	Corn
End mass	0.506 g
Humidity content	14.237% M

Example of a complex drying report, generated by PMV 50.5Y moisture analyzer.

# Databases Managing and Editing

## Databases Drying Process Ergonomics

Drying parameters such as temperature and automatic shutdown are optimally selected for every single sample. Trying to remember the parameters for just a few samples requires considerable effort. It is more convenient to record the parameters in a database than attempting to learn them by heart.

### Product Database Contains any Data Relating to a Sample:

- name and description,
- EAN code: searching a sample in a database by means of a scanner,
- target value (%): value used for automatic control of sample weight (bar graph) and for determining moisture content variation over time (trend),
- Min, Max: value used for automatic control of sample weight (bar graph),
- tolerance: value used for determining moisture content variation over time (trend),
- drying program.

### Drying Programs Database Contains any Data Relating to a Drying Process:

- name, code,
- drying mode, drying temperature,
- automatic shutdown (auto / time-defined / user-defined),
- start mass control (none / optional / essential),
- equipment intended for a sample (methodology),
- instruction on sample preparation for a drying process (methodology),
- required sample size (methodology).

## Database Editor


Database Editor PC software is designed to support users dealing with a vast number of samples. Clear structure of the program ensures quickness when it comes to specifying drying parameters and other information relating to a sample. Data is transferred from the software to a moisture analyzer by means of Ethernet ( ) or RS 232 ( , R).

Users	Code	Name	Desc
	12	Mustard	Musta
	13	Powdered Soya Drink	Powd
	14	Pistachio Nut	Pistac
	15	Walnut	Waln
	16	Wheat Bran	Whea
	17	PA 6	PA 6
	18	Fodder	Fodd
	19	Soya Pate With Mushrooms	Soya
	20	PC (polycarbonate)	PC (p
	21	Pellet	Pellet
	22	Gingerbread	Ginge
	23	Corn Flakes	Corn
	24	Washing-up Detergent	Wash

Available databases: Products, Weighing Records, Customers, Drying Programs, Drying Process Records, Ambient Conditions, Packaging, Warehouses, Printouts, Universal Variables.

Name:	Walnut		
Description:	Walnut		
Code:	3	Code EAN:	9854327
Target value:	4,23	Unit:	%M
Drying program:		Tare:	0
Min:		Max:	4,5
Tolerance:	Powdered Soya Drink Pistachio Nut <b>Walnut</b> Wheat Bran Candied Papaya Granulated Fodder Soya Pate With Mushrooms		
Density:		Shelf life-time:	
Price:		VAT:	
Date:			
Printout:	✕		

Detailed information concerning the product.




Export / import of databases moisture analyzers

# E2R Results Analysis

## Databases Drying Processes Ergonomics

Production processes for which moisture content of a particular sample is a crucial parameter, require quick reaction. This may be carried out using so called weighing networks comprising PMV 50 PLUS, MA X2 and MA R moisture analyzers. Each drying process is monitored on-line regardless of workstations location.

E2R Moisture Analyzer PC software is designed to record measurements performed by means of RADWAG moisture analyzers cooperating in a network, using RS 232 and Ethernet interfaces for connection. the software enables monitoring and reporting of collected measurements.



Software functions: on-line monitoring of moisture analyzer operation, possibility of configuring reports and graphs, analysis of data collected from many drying workstations, data protection.

## PC Software

RADWAG PC software supports moisture analyzers expanding their functionality.

### R-Lab

Scales preview, weighings graphs and statistics graphs.

### Database Editor

Readout, databases editing and record of computer stored databases on balance.

### Rad Key

Readout of balance data by means of defined Hot Key.

### E2R Moisture Analyzers

Record of weighments carried out by moisture analyzers cooperating in a network.

## Additional Equipment

- Anti-vibration weighing tables,
- Disposable weighing pans,
- Thermal and dot matrix printers,
- Barcode scanners,
- Control thermometer,
- Water vapor permeability set.

Complete offer is to be found on [www.radwag.com](http://www.radwag.com) website.



### E2R Moisture Analyzer is

**a module of integrated system for managing E2R weighing processes.** E2R System comprises various programs ensuring continuous control of balances and their databases together with both, complete managing of the manufacturing process and the process optimisation.



# Intended Use and Applications

## Area of Use

Moisture content analysis and dry mass measurement of a particular product are both crucial for various branches of industry and science. Vast area of use and diversity of analyzed samples structure require individual approach to different substances.



### Dairy Industry

**Samples:** cheese, buttermilk, yoghurt, powdered milk, etc. Samples are dried directly on a weighing pan or by means of glass fiber filters or silica sand (increasing surface of evaporation).



### Fruit and Vegetable Industry

**Samples:** dried vegetables, fruits and mushroom, nuts etc. Samples shall be cut into smaller pieces (the analyzed samples cannot be too thick).



### Food Industry

**Samples:** sugar, flour, pasta, spices, gelatin, etc. Thin layer of semi-liquid samples shall be distributed on a weighing pan (silica sand or glass fiber filters may be used). Other kinds of samples shall be crushed.



### Chemical Industry

**Samples:** emulsion, gel and lotions used for cleaning, paints, film, graphite, etc. Thin layer of semi-liquid samples shall be distributed on a weighing pan (Silica sand or glass fiber filters may be used). Other kinds of samples shall be crushed.



### Agricultural Industry

**Samples:** grain, seeds, hay, biomass, etc. Grain needs to be crushed prior drying.

# Technical Specification

## Microwave Moisture Analyzer



Product Code	Model	Max. capacity	Readability	Heating module
WL-307-0006	<b>PMV 50.5Y</b>	50 g	0,1 mg	microwave radiation emitter

## X2 Series



Product Code	Model	Max. capacity	Readability	Heating module
WL-306-0032	<b>MA 50.X2.A.WH</b>	50 g	1 mg	halogen
WL-306-0043	<b>MA 50.X2.IC.A.WH</b>	50 g	1 mg	halogen
WL-306-0035	<b>MA 50/1.X2.A.WH</b>	50 g	0,1 mg	halogen
WL-306-0047	<b>MA 50/1.X2.IC.A.WH</b>	50 g	0,1 mg	halogen
WL-306-0021	<b>MA 50.X2.A</b>	50 g	1 mg	IR emitter
WL-306-0041	<b>MA 50.X2.IC.A</b>	50 g	1 mg	IR emitter
WL-306-0026	<b>MA 50/1.X2.A</b>	50 g	0,1 mg	IR emitter
WL-306-0045	<b>MA 50/1.X2.IC.A</b>	50 g	0,1 mg	IR emitter
WL-306-0063	<b>MA 50.X2.A.NS</b>	50 g	1 mg	metal heater
WL-306-0074	<b>MA 50.X2.IC.A.NS</b>	50 g	1 mg	metal heater
WL-306-0075	<b>MA 50/1.X2.A.NS</b>	50 g	0,1 mg	metal heater
WL-306-0076	<b>MA 50/1.X2.IC.A.NS</b>	50 g	0,1 mg	metal heater
WL-306-0034	<b>MA 110.X2.A.WH</b>	110 g	1 mg	halogen
WL-306-0051	<b>MA 110.X2.IC.A.WH</b>	110 g	1 mg	halogen
WL-306-0027	<b>MA 110.X2.A</b>	110 g	1 mg	IR emitter
WL-306-0049	<b>MA 110.X2.IC.A</b>	110 g	1 mg	IR emitter
WL-306-0080	<b>MA 110.X2.A.NS</b>	110 g	1 mg	metal heater
WL-306-0081	<b>MA 110.X2.IC.A.NS</b>	110 g	1 mg	metal heater
WL-306-0078	<b>MA 200/1.X2.A.WH</b>	200 g	0,1 mg	halogen
WL-306-0065	<b>MA 200/1.X2.IC.A.WH</b>	200 g	0,1 mg	halogen
WL-306-0061	<b>MA 200/1.X2.A</b>	200 g	0,1 mg	IR emitter
WL-306-0062	<b>MA 200/1.X2.IC.A</b>	200 g	0,1 mg	IR emitter
WL-306-0077	<b>MA 200/1.X2.A.NS</b>	200 g	0,1 mg	metal heater
WL-306-0079	<b>MA 200/1.X2.IC.A.NS</b>	200 g	0,1 mg	metal heater
WL-306-0033	<b>MA 210.X2.A.WH</b>	210 g	1 mg	halogen
WL-306-0055	<b>MA 210.X2.IC.A.WH</b>	210 g	1 mg	halogen
WL-306-0028	<b>MA 210.X2.A</b>	210 g	1 mg	IR emitter
WL-306-0053	<b>MA 210.X2.IC.A</b>	210 g	1 mg	IR emitter
WL-306-0082	<b>MA 210.X2.A.NS</b>	210 g	1 mg	metal heater
WL-306-0083	<b>MA 210.X2.IC.A.NS</b>	210 g	1 mg	metal heater

## R Series



Product Code	Model	Max. capacity	Readability	Heating module
WL-305-0005	<b>MA 50.R.WH</b>	50 g	1 mg	halogen
WL-305-0006	<b>MA 50/1.R.WH</b>	50 g	0,1 mg	halogen
WL-305-0001	<b>MA 50.R</b>	50 g	1 mg	IR emitter
WL-305-0002	<b>MA 50/1.R</b>	50 g	0,1 mg	IR emitter
WL-305-0014	<b>MA 50.R.NS</b>	50 g	1 mg	metal heater
WL-305-0039	<b>MA 50/1.R.NS</b>	50 g	0,1 mg	metal heater
WL-305-0007	<b>MA 110.R.WH</b>	110 g	1 mg	halogen
WL-305-0004	<b>MA 110.R</b>	110 g	1 mg	IR emitter
WL-305-0024	<b>MA 110.R.NS</b>	110 g	1 mg	metal heater
WL-305-0008	<b>MA 210.R.WH</b>	210 g	1 mg	halogen
WL-305-0003	<b>MA 210.R</b>	210 g	1 mg	IR emitter
WL-305-0009	<b>MA 210.R.NS</b>	210 g	1 mg	metal heater
WL-305-0009	<b>MA 210.R.NS</b>	210 g	1 mg	metal heater