



PRODUCT CATALOGUE 2022



radwag

Weighing Precision John Smith Wed, 10.02.2022 09:31:15

0.00000 g DWA


-I- 0 Max

Supplementary unit

Gross	0.00000 g	Temperature	19°C	Air Humidity	45%
Tare	0.00000 g	Vibrations	3%	Air Pressure	980.1 hPa
Product	Liquid	Air Density	1.158 kg/m ³		

radwag.com

Adherence to law regulations
Information safety
Intuitive operation




www.radwag.com

PGC


Prepackaged Goods Control

STATISTICAL CONTROL


Law-accordant and internal control criteria Non-destructive and destructive control




PGC control, HY10.PGC.H



HY10.PGC.H Scale's Home Screen



PGC weighings carried out by means of HY10.PGC.H scale



Control completion screen on WLY/PGC and WPY/PGC scales

Versatile Use of PGC Control Scales

PGC control scales are used in food, chemical, pharmaceutical and cosmetic industry wherever statistical control is a must.

Adherence to International Law Regulations

Software accordant with Polish law on prepackaged goods control, PGC Law of 7th May 2009 (OJ No. 91, item 740), Directive 76/211/EEC and WELMEC 6.4 Guide titled "Guide for Packers and Importers of e-marked Prepacked Products".

Innovative Solutions for Control Scales

PGC control scales have been equipped with high-tech colour touchscreen providing perfect readability even in low light conditions. Clear menu arrangement improves performance and adds to maximum comfort of operation.

Independent Operation of PGC Scales

PGC scales operation is based on local database whereas computer database is a source of information on operators, products and time schedules for scales cooperating in a network.

Authorized Access Control

Access levels to particular functions, defined individually for particular users, prevent interference of unauthorized operators.

Databases Management

Option of adding, deleting, and editing local products, operators and time schedules database.

Safe Mechanisms for Data Record and Storage

Upon control completion final report featuring all necessary data is saved. Control and average tares reports can be printed or recorded to USB flash drive.

User-friendly Interface

At each control stage important control-related information is displayed, with this comfort of operation is guaranteed.

Optimization and Production Cost Reduction

Declaring tolerance high threshold allowing to maintain control over product giveaway.

Internal Criteria

PGC scales offer option of carrying out control in accordance with internal criteria, including tolerance thresholds. The said option enables adaptation of the whole control process to internal regulations and operations performed by the customer.

Various Methods of Batch Assessment

Using the scale it is possible to perform control by means of few methods (non-destructive average tare, non-destructive empty-full and destructive full-empty).

Performance Boost

The scale facilitates performance of two concurrent controls, restarted upon power supply loss.

Controls Schedule

Going according to a predefined control schedule or cyclic control of selected products within specified time interval.

Auto Reminders

By means of prompts displayed on a screen or on an external signaling device the scale reminds you about necessity to carry out a measurement in the course of control. You are also prompted when it is time to start the control.

Average Tare Estimation Module

An in-built tare estimation module enables performance of average tare control directly prior product control, it also facilitates cyclic supervision within predefined time interval.

PC Software PGC Viewer




Diagram of data acquisition carried out by means of PGC Viewer

PGC Viewer is an application designed to support workstations for prepackaged goods control. Using the application you can preview and print reports on accomplished controls of prepackaged goods

Defined column	
Controls	
	CONTROL REPORT NO. U/07/07/06/13/16
	CONTROL REPORT NO. U/07/07/06/13/16
	CONTROL REPORT NO. Z/07/07/06/15/55

PC Software PGC Viewer

E2R PGC An integrated module




Diagram of E2R PGC system


Modular E2R Software

E2R PGC module of a modern, fully integrated weighing system, E2R, is a tool providing you with maximum security. The system, due to innovative modules characteristic for their extensive functionality, enables automation, full support and control of the manufacturing processes.


Advanced functions facilitate monitoring of current devices state, remote control start and preview of complex analysis, i.e. reports and charts, of archive controls.

Cutting-edge technology is a warranty of faultless system operation and stored data safety. Data recorded on a balance is sent to a computer database based on Microsoft SQL server. By means of E2R PGC application it is possible to control particular users access to the system.

When it comes to exchange of data, the system offers option of cooperation with external IT systems. All the above plus user-friendly environment add to your comfort.



PUE 7.1 Home Screen



Online preview of measurements



Scales for Statistical Control - Wide Range of Maximum Capacities
Max capacity of scales for statistical control ranges from 0.6g to 120kg, accuracy starts from 10mg.

Reliability and Hygiene in Challenging Working Conditions
Stainless steel mechanical design and IP68/69 in-use facilitate operation of HY10.PGC.H under challenging industrial conditions such as great dust and high humidity. These two characteristics make the scales meet high hygiene standards for food and pharmaceutical industries.

Ergonomics and Comfort of Operation
Due to compact and ergonomic design, access to the interfaces of WPY/PGC scales is convenient. The device can be operated using barcode reader, transponder card reader and IR sensors. This enables touch free operation of the workstation.

Redefined Accuracy
WLY/PGC series is intended for fast and precise mass measurement under both laboratory and industrial conditions, wherein the readability is even 0.01g.

Professional Weighing and the Highest Measurement Accuracy
With the highest possible weighing accuracy and measurement repeatability of $d < 1$, redefines quality of mass measurement. Control over the right scale's level is carried out by means of a semi-automatic levelling system. An in-built system for ambient conditions monitoring facilitates control of the workroom. As for the measurement accuracy it is provided with auto adjustment function.

Technical Specification

Packaged Goods Control Scales - Statistic Control

PS X2 Balances for PGC



Product Code	Model	Max. capacity	Readability	Weighing pan dimensions	Verification scale interval
WK-306-0005	PS 750.X2.PGC	750 g	0.01 g	128x128 mm	0,01 g



Product Code	Model	Max. capacity	Readability	Weighing pan dimensions	Verification scale interval
WK-306-0006	PS 2100.X2.M.PGC	2100 g	0.1 g	195x195 mm	0,1 g
WK-306-0007	PS 4500.X2.M.PGC	4500 g	0.1 g	195x195 mm	0,1 g
WK-306-0008	PS 6100.X2.M.PGC	6100 g	0.1 g	195x195 mm	0,1 g
WK-306-0009	PS 8100.X2.M.PGC	8100 g	0.1 g	195x195 mm	0,1 g

WLY Scales for PGC



Product Code	Model	Max. capacity	Readability	Weighing pan dimensions	Verification scale interval
WK-307-1034	WLY 0.6/1.2/PGC/D2	0.6 / 1.2 kg	0.1 / 0.2 g	195×195 mm	0,1 / 0,2 g
WK-307-1001	WLY 3/6/PGC/D2	3 / 6 kg	0.5 / 1 g	195×195 mm	0,5 / 1 g
WK-307-1040	WLY 6/12/PGC/D2	6 / 12 kg	1 / 2 g	195×195 mm	1 / 2 g



Product Code	Model	Max. capacity	Readability	Weighing pan dimensions	Verification scale interval
WK-307-1042	WLY 6/12/PGC/F1/R	6 / 12 kg	1 / 2 g	300×300 mm	1 / 2 g
WK-307-1041	WLY 6/12/PGC/F1/K	6 / 12 kg	1 / 2 g	300×300 mm	1 / 2 g



Product Code	Model	Max. capacity	Readability	Weighing pan dimensions	Verification scale interval
WK-307-0002	WLY 60/PGC/C2/K	60 kg	10 g	400×500 mm	10 g
WK-307-0001	WLY 60/PGC/C2/R	60 kg	10 g	400×500 mm	10 g
WK-307-1022	WLY 60/120/PGC/C2/K	60 / 120 kg	10 / 20 g	400×500 mm	10 / 20 g
WK-307-1021	WLY 60/120/PGC/C2/R	60 / 120 kg	10 / 20 g	400×500 mm	10 / 20 g

WPY Scales for PGC



Product Code	Model	Max. capacity	Readability	Weighing pan dimensions	Verification scale interval
WK-308-1046	WPY 0.6/1.5/PGC/D2	0.6 / 1.5 kg	0.2 / 0.5 g	195×195 mm	0.2 / 0.5 g
WK-308-1002	WPY 1,5/3/PGC/D2	1.5 / 3 kg	0.5 / 1 g	195×195 mm	0.5 / 1 g
WK-308-1003	WPY 3/6/PGC/D2	3 / 6 kg	1 / 2 g	195×195 mm	1 / 2 g



Product Code	Model	Max. capacity	Readability	Weighing pan dimensions	Verification scale interval
WK-308-1049	WPY 3/6/PGC/F1/K	3 / 6 kg	1 / 2 g	300×300 mm	1 / 2 g
WK-308-1050	WPY 3/6/PGC/F1/R	3 / 6 kg	1 / 2 g	300×300 mm	1 / 2 g
WK-308-1051	WPY 6/15/PGC/F1/K	6 / 15 kg	2 / 5 g	300×300 mm	2 / 5 g
WK-308-1052	WPY 6/15/PGC/F1/R	6 / 15 kg	2 / 5 g	300×300 mm	2 / 5 g
WK-308-1053	WPY 15/30/PGC/F1/K	15 / 30 kg	5 / 10 g	300×300 mm	5 / 10 g
WK-308-1054	WPY 15/30/PGC/F1/R	15 / 30 kg	5 / 10 g	300×300 mm	5 / 10 g



Product Code	Model	Max. capacity	Readability	Weighing pan dimensions	Verification scale interval
WK-308-1035	WPY 15/30/PGC/C2/K	15 / 30 kg	5 / 10 g	400×500 mm	5 / 10 g
WK-308-1036	WPY 15/30/PGC/C2/R	15 / 30 kg	5 / 10 g	400×500 mm	5 / 10 g
WK-308-1034	WPY 30/60/PGC/C2/K	30 / 60 kg	10 / 20 g	400×500 mm	10 / 20 g
WK-308-1033	WPY 30/60/PGC/C2/R	30 / 60 kg	10 / 20 g	400×500 mm	10 / 20 g

HY10.H Scales for PGC




Product Code	Model	Max. capacity	Readability	Weighing pan dimensions	Verification scale interval
WK-311-1002	HY10.1,5/3.PGC.H1	1.5 / 3 kg	0.5 / 1 g	150×200 mm	0,5 / 1 g
WK-311-1004	HY10.3/6.PGC.H1	3 / 6 kg	1 / 2 g	150×200 mm	1 / 2 g
WK-311-1006	HY10.3/6.PGC.H2	3 / 6 kg	1 / 2 g	250×300 mm	1 / 2 g
WK-311-1007	HY10.6/15.PGC.H2	6 / 15 kg	2 / 5 g	250×300 mm	2 / 5 g
WK-311-1008	HY10.6/15.PGC.H3	6 / 15 kg	2 / 5 g	410×410 mm	2 / 5 g
WK-311-1009	HY10.15/30.PGC.H3	15 / 30 kg	5 / 10 g	410×410 mm	5 / 10 g
WK-311-1005	HY10.30/60.PGC.H3	30 / 60 kg	10 / 20 g	410×410 mm	10 / 20 g



PRODUCT CATALOGUE

ADVANCED WEIGHING TECHNOLOGIES





www.radwag.com