

GAS ANALYZER

STATIONARY GAS ANALYSERS AND GAS ANALYSIS SYSTEMS FOR LABORATORY AND INDUSTRY



- Multi-channel systems for higher economic efficiency
- Infrared, paramagnetic, electrochemical and PID gas sensors:
SO₂, O₂, CO₂, CO, CH₄, H₂S, NH₃, SO₂, VOC
- High precision sensors – from 100% to ppm levels
- Serial interface USB/RS232, RS485, ETHERNET
- Touch screen interface
- Auto-calibration option
- 1 to 8 sensors in 19” enclosure



Accredited laboratory for gas flow measurements SIST EN ISO 17025:2005.

MANUFACTURER

Echo d.o.o.
Stari trg 37, Slovenske Konjice 3210

tel: +386 (0)3 759 23 80
info@echo.si
www.echoinstruments.eu

DISTRIBUTOR

