



PRODUCT CATALOGUE 2022





www.radwag.com

Software

ADVANCED WEIGHING TECHNOLOGIES

E2R System



Collect data from moisture analyzers that are accessible via the company network. You can have up to 16 devices under constant control. Analyze results, make graphs and reports, archive data.

Read more:



R-Cloud




The R-Cloud is a RADWAG service that enables access to results stored in a cloud. Access to the stored data is possible from any place with Internet connection regardless of the location. Open <https://rcloud.radwag.com/> website and carry out registration procedure.

Read more:



R-Lab




The modern R-Lab software enables collecting, presenting and carrying out statistical analysis of measurements sent from RADWAG scales to a computer. the software operates with max 20 scales simultaneously.

Read more:



RMCS Lite



RMCS Lite software is designed to operate RADWAG-manufactured mass comparators via a computer working in one and the same local network. the software allows for carrying out mass measurement in accordance with PN-EN 12341:2014 standard.

Read more:



and more...

Do you want to see all our software?

Visit our website www.radwag.com, or use this QR Code.





Contact Us

Sales

Laboratory Balances

+48 483 866 314
m.staniak@radwag.com

Technical support

labbalances@radwag.com

Scales for Ex Areas

Sales
+48 483 866 310
g.rog@radwag.pl

Technical support

atex@radwag.com

Mass Comparators

Sales
+48 483 866 314
m.staniak@radwag.com

Technical support

masscomparators@radwag.com

Others

Sales
+48 483 866 314
m.staniak@radwag.com

Technical support

validation@radwag.com

Industrial Scales

Sales
+48 483 866 310
g.rog@radwag.pl

Technical support

indscales@radwag.com

Checkweighers

Sales
+48 483 866 310
g.rog@radwag.pl

Technical support


checkweighers@radwag.com

Software

Sales
+48 483 866 310
g.rog@radwag.pl

Technical support

pcsoftware@radwag.com



www.radwag.com